



## eCTAS Release Notes: Version 2.33.0.0

Release Date: Dec 5, 2024

Integrations: **B** - Basic **C** - Complex **WS** - Web Service **CRT** - Certification

Change Type	Details	B	C	WS	CRT
User Interface (UI) Changes	<p><b>Triage in Progress indicator appears immediately</b></p> <ul style="list-style-type: none"> <li>When a new patient is received in eCTAS and the nurse selects the 'Continue Triage' button, the brief delay between the patient's name appearing on the Pretriaged Patients queue and the Triage in Progress indicator being displayed has been removed. The Triage in Progress indicator is now visible as soon as the queue refreshes.</li> </ul>	✓	✓		
	<p><b>Test Result Report Tool</b></p> <ul style="list-style-type: none"> <li>A new option has been added that allows users to send an automated email to Ontario Health when they have results and / or comments ready for review.</li> </ul>				✓
	<p><b>Administration Console</b></p> <ul style="list-style-type: none"> <li>A tool that allows Ontario Health administrators to manage nurse notifications has been enhanced. This tool is not visible to hospital users.</li> </ul>	N/A	N/A	N/A	N/A
HL7 Message Processing	<p><b>Trailing spaces in modifier names</b></p> <ul style="list-style-type: none"> <li>Trailing spaces that appeared at the end of some modifier names in HL7 ORU Triage messages have been removed.</li> </ul>		✓		

Change Type	Details	B	C	WS	CRT
	<p><b>CTAS Guideline Change Readiness</b></p> <ul style="list-style-type: none"> <li>In preparation for the upcoming CTAS Guideline changes, minor updates have been made to HL7 message processing to ensure that the current and upcoming versions of the guidelines can co-exist during the transition period.</li> <li>These updates <u>do not</u> impact the structure or content of any existing HL7 messages.</li> </ul>		✓		
Web Service Enhancements	<p><b>SpO2 Method Validation</b></p> <ul style="list-style-type: none"> <li>The validation error indicating that the submitted SpO2 method value is invalid has been changed to a warning. Triage containing an invalid value in this field will now be saved with the invalid value excluded.</li> </ul>			✓	✓
Bug fixes	<p><b>eCTAS Clinical Application</b></p> <ul style="list-style-type: none"> <li>A workflow issue that prevented the Triage in Progress indicator from being displayed in rare cases has been corrected.</li> <li>A synchronization issue that resulted in Priority in Queue notes and TOA details not refreshing in the Summary panel next to the patient queues has been corrected.</li> <li>The triage print report now lists Infection Control Screening (ICS) responses in the same order that the questions are displayed on the ICS screen regardless of the order answers are entered during screening.</li> <li>A configuration issue that prevented some notification messages from displaying at their scheduled time has been corrected.</li> </ul>	✓	✓		
	<p><b>HL7 message processing fixes</b></p> <ul style="list-style-type: none"> <li>An issue that caused rare message processing failures due to a missing message acknowledgement value has been resolved.</li> </ul>		✓		

Change Type	Details	B	C	WS	CRT
	<p><b>Administration Console</b></p> <ul style="list-style-type: none"> <li>• A configuration issue that prevented the 'Sign Out' link from returning users back to the Sign In page has been corrected.</li> <li>• Incorrect data shown in the Site Role Management tool when a user was granted access to 2 or more new sites but only assigned roles at one of the sites has been corrected.</li> </ul>	✓	✓	✓	✓
Performance Improvements	<p><b>System upgrades</b></p> <ul style="list-style-type: none"> <li>• System component upgrades and instrumentation enhancements have been made to ensure ongoing system stability and security.</li> </ul>	✓	✓	✓	✓