

## **Release Notes: Version 2.6.0.0**

Release Date: May 22, 2018



Note: Nurses may need to <u>clear their browser cache</u> if changes are not immediately visible in the eCTAS application.

Integrations: B - Basic C - Complex WS - Web Service CRT - Certification

Change Type	Details	В	С	ws	CRT
Performance Improvement	Users will notice an improvement associated with syncing of patient information between workstations. Users should experience a greater level of consistency in their overall experience.	<b>⊘</b>	<b>⊘</b>		
	<ul> <li>HL7 Messages</li> <li>For messages transmitted to eCTAS from a Hospital Information System, users should notice increased performance associated with displaying this information within the eCTAS application.</li> </ul>				
	The Refresh Button  The performance of the Refresh button has been improved to increase the level of consistency. The frequency in which a user is required to use the refresh function should be greatly reduced.  Refresh button				
Bug Fix	Changes to patient queues were previously not being updated upon logging in to eCTAS for some users. This issue has been resolved, the Patient Queues will display an up-to-date list of patient records upon logging in for all users.	<b>Ø</b>	<b>⊘</b>		
	HL7 message failures caused by invalid characters associated to vital signs have been resolved. The system will prevent the copy and paste of invalid characters and will only accept numerical values within vital sign fields. This will reduce instances of failures associated with the transmission of triage and reassessment messages from eCTAS to your Hospital Information Systems.				