

eCTAS Release Notes: Version 2.37.0.0

Release Date: Apr 7, 2026

Integrations: **B** - Basic **C** - Complex **WS** - Web Service **CRT** - Certification

Change Type	Details	B	C	WS	CRT
Hospital Requested Enhancement	Health card barcode scanning <ul style="list-style-type: none"> Nurses receiving and logging patients into eCTAS can now use a barcode scanner to read data from a patient's health card as an alternative to a magnetic stripe reader. Hospitals are responsible for acquiring barcode scanners to make use of this feature. Scanners must support PDF417 format. 	✓	✓		
	Triage record printout update <ul style="list-style-type: none"> The layout of the triage record printout has been updated to remove the large blank spaces that could appear between sections of the report. This improves readability and limits the number of pages required to print out the triage record. The layout changes apply for triages printed from both the eCTAS Clinical Application and the Registration View. 	✓	✓	✓	✓
User Interface (UI) Changes	Visual indicator for automatically selected CTAS 1 scores <ul style="list-style-type: none"> A new visual indicator has been added on the Triage Assessment screen to alert nurses when an automatically selected modifier has updated the triage score to a CTAS 1. The purpose of the indicator is to encourage nurses to review whether CTAS 1 is a clinically appropriate score for the patient and manually adjust modifiers as necessary. Nurses are not required to take any action to acknowledge or dismiss this indicator. It is purely a visual cue to review the applied CTAS score. 	✓	✓		

Change Type	Details	B	C	WS	CRT
	<p>Test Result Report Tool</p> <ul style="list-style-type: none"> The Test Result Report tool has been expanded to allow Application Complex hospitals to view the status of their Conformance Test results. Facility Role Administrators can assign the 'Site Integration Tester' role to any testers who need access to the tool. The role name has changed from 'Site API Tester'. All users previously assigned access will be automatically updated to the new role name. The Test Result Report tool is only available in the Conformance environment. 		✓		
	<p>Administration Console</p> <ul style="list-style-type: none"> The tool that allows Ontario Health administrators to manage Conformance and Clinical Validation test phases has been updated. This tool is not visible to hospital users. 	N/A	N/A	N/A	N/A
Bug fixes	<p>Triage Printout</p> <ul style="list-style-type: none"> A few minor formatting issues with the triage printout have been corrected: <ul style="list-style-type: none"> A blank line that appeared between the SpO2 value and the selected method when no source was selected on a reassessment has been removed. An empty EMS notes section that appeared when the mode of arrival was not EMS has been removed. Units are now displayed if a CBGM value of 0 is entered. Units are now displayed if a weight value of 0 is entered. 	✓	✓	✓	✓
	<p>Triaged Patients Queue</p> <ul style="list-style-type: none"> A rare display error that could result in the 'Saving Assessment' animation to not disappear a few seconds after clicking the 'I'm Done Triage' button has been corrected. 	✓	✓		

Change Type	Details	B	C	WS	CRT
	<p>Previous Visit Lookup</p> <ul style="list-style-type: none"> The previous visit lookup has been optimized to prevent lookup attempts when the patient date of birth is not provided. 	✓	✓		
Performance Improvements	<p>System upgrades</p> <ul style="list-style-type: none"> System component upgrades and instrumentation enhancements have been made to ensure ongoing system stability and security. 	✓	✓	✓	✓